



PART 1: Installation of extension « MRI_Quality_Control » in 3DSlicer

- **To do only once**

- Launch 3dSlicer
- 1) « Welcome to Slicer » and « All Modules »
- 2) « Extension Wizard »



- 3) « Select Extension »




- Select the file « MRI_Quality_Control » available here (need to be unzip before):

https://github.com/AurelienCD/MRI_Quality_Control/blob/master/MRI_Quality_Control.zip

- Restart 3DSlicer and then, extension will be available all the time and will not need to be loaded as before

PART 2: How to use « MRI_Quality_Control » extension ?

- Load CT image of the phantom with the different labels (available on request: a.corroyer-dulmont@baclesse.unicancer.fr)
=> link: <https://owncloud.cyceron.fr/index.php/s/N5xRiMA61zoptNk/download>
- Load the MRI image
- 1) « Welcome to Slicer »
- 2) « Imaging »
- 3) « MRI_Quality_Control »




- Reference CT and spheres labels are normally automatically affected, just check it and then chose the good image for the MRI

Parameters	
Reference CT Image:	CT_REF
Spheres label (CT space):	Spheres_label
MR Image to analyse:	A1.bac.be6b.fr.badesse.S1.5993684.1_00000.DCM

- Chose the directory output where the results will be saved
- Follow the different steps as in the picture below. Wait the end of each step with the green message at the bottom of the soft in the python interactor window.

1) Registration
2A) Geometric distortion analysis
Geometric distortion result image (only for testing, let it in none):
2B) Show 3D Geometric distortion results
3) Signal to Noise Ratio analysis
4) Signal Linearity analysis
5) Geometric Accuracy analysis

- Click on « **1) Registration** » to launch the registration between the CT ref and the MRI
- When registration is finished, click on « **2A) Geometric distortion analysis** »
- When the analyse is done, click on « **2B) Show 3D Geometric distortion results** » to obtain the 3D representation of the distortion in the spheres. Colors represent the geometric distortion (green<1mm; orange<2mm; red>2mm)



- « **3) Signal to Noise Ratio analysis** »
- « **4) Signal Linearity analysis** »
- « **5) Geometric Accuracy analysis** »
- All the results are in the output directory in a text file